


Assut sutures

of Switzerland


Vascular & Cardiovascular Sutures



A leading Swiss brand, well-known worldwide

Astralen/Polyester coated

The sterile Assut Astralen coated suture consists of braided polyester with a fine and uniform coating of medical grade silicone. A modern braiding technique and a special coating process provide a soft and easy to handle suture material with high tensile strength and knot security. Non-absorbable and available in green and white, single and multipack version; a special assortment is pre-threaded with pledgets.



Astralen (polyester)
braided coated



Needle code	USP	EP	CM	Product No. (REF)
½ CIRCLE - TAPER POINT (ROUND-BODIED)				
2 x HR-17 	2-0	3	75cm	CV 3875-35x10 5 green & 5 white braided sutures with 3x5x1.5mm oval firm pledget
	2-0	3	75cm	CV 3875-33x10 5 green & 5 white braided sutures with 3x3x1.5mm oval firm pledget
	2-0	3	75cm	CV 3875-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
	2-0	3	75cm	CV 3875x10 5 green & 5 white braided sutures
	2-0	3	75cm	CV 38855-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
	2-0	3	75cm	CV 38855x10 5 green & 5 white braided sutures
	2-0	3	75cm	CV 3945-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
2 x HR-25 	2-0	3	75cm	CV 3945-24x6 3 green & 3 white braided sutures with 2x4x1.5mm rectangular firm pledget
	2-0	3	75cm	CV 3945x10 5 green & 5 white braided sutures
	2-0	3	75cm	CV 39485-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
	2-0	3	75cm	CV 39485x10 5 green & 5 white braided sutures
	2-0	3	75cm	CV 39485x10 5 green & 5 white braided sutures
½ CIRCLE - CUTTING TAPER				
2 x HSR-18 	2-0	3	75cm	CV 3945-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
	2-0	3	75cm	CV 3945-24x6 3 green & 3 white braided sutures with 2x4x1.5mm rectangular firm pledget
	2-0	3	75cm	CV 3945x10 5 green & 5 white braided sutures
	2-0	3	75cm	CV 39485-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
2 x HSR-25 	2-0	3	75cm	CV 39485-37x10 5 green & 5 white braided sutures with 3x7x1.5mm oval firm pledget
	2-0	3	75cm	CV 39485x10 5 green & 5 white braided sutures
	2-0	3	75cm	CV 39485x10 5 green & 5 white braided sutures

Astralen/Polyester coated







Needle code	USP	EP	CM	Product No. (REF)
½ CIRCLE - TAPER POINT (ROUND-BODIED)				
2 x HR-17 	4 - 0	1.5	75 cm (green)	901
	4 - 0	1.5	90 cm (green)	901N
	3 - 0	2	75 cm (green)	902
	3 - 0	2	90 cm (green)	902N
	2 - 0	3	75 cm (green)	3875
	2 - 0	3	90 cm (green)	3875N
2 x HR-25 	3 - 0	2	75 cm (green)	38854
	3 - 0	2	90 cm (green)	38854N
	2 - 0	3	75 cm (green)	38855
	2 - 0	3	90 cm (green)	38855N
2 x HR-27 	3 - 0	2	75 cm (green)	903
	3 - 0	3	90 cm (green)	903N
	2 - 0	3	75 cm (green)	904
	2 - 0	3	90 cm (green)	904N
2 x HR-35 	2 - 0	3	90 cm (green)	38975N
½ CIRCLE - CUTTING TAPER				
2 x HSR-18 	5 - 0	1	90 cm (green)	3942N
	4 - 0	1.5	90 cm (green)	3943N
	3 - 0	2	90 cm (green)	3944N
	2 - 0	3	90 cm (green)	3945N
	2 - 0	3	90 cm (green)	3945Nx4
	Multipack 4 sutures per envelope			
2 x HSR-25 	2 - 0	3	75 cm (green)	39485
	2 - 0	3	90 cm (green)	39485N
2 x HSR-27 	3 - 0	2	90 cm (green)	3954N
	2 - 0	3	90 cm (green)	3955N
	- 0 -	3.5	90 cm (green)	3956N
	2 - 0	3	90 cm (green)	3955Nx4
	Multipack 4 sutures per envelope			

Polypropylene blue


A blue monofilament sterile suture material with exceptional smooth and uniform surface allowing easy passage through tissue with minimum trauma. Extremely well tolerated, this non-absorbable suture has excellent knot security.



Polypropylene blue monofilament

Needle code	USP	EP	CM	Microns	Description	Product No. (REF)
3/8 CIRCLE - TAPER POINT (ROUND-BODIED)						
2 x TRM-6,5 	8 - 0	0.4	60 cm	140 Mi.	6,5 mm micropoint, double-armed	91051408
2 x TRM-8 	8 - 0	0.4	60 cm	150Mi.	8 mm micropoint, double-armed	950508
	7 - 0	0.5	60 cm	150 Mi.	8 mm micropoint, double-armed	950507
2 x TRM-9,3 	7 - 0	0.5	60 cm	200 Mi.	9,3 mm micropoint, double-armed	952607
	7 - 0	0.5	75 cm	200 Mi.	9,3 mm micropoint, double-armed	952697
2 x TRM-10 	7 - 0	0.5	60 cm	200 Mi.	10 mm micropoint, double-armed	954607
	7 - 0	0.5	75 cm	200 Mi.	10 mm micropoint, double-armed	954697
2 x TRM-13 	7 - 0	0.5	60 cm	230 Mi.	13 mm micropoint, double-armed	960707
1/2 CIRCLE - TAPER POINT (ROUND-BODIED)						
2 x HRM-10 	7 - 0	0.5	75 cm	280 Mi.	10 mm micropoint, double-armed	910697
3/8 CIRCLE - CUTTING TAPER						
2 x TSRM-9 	7 - 0	0.5	75 cm	230 Mi.	9 mm micropoint, double-armed	972797
2 x TSRM-13 	6 - 0	0.7	75 cm	300 Mi.	13 mm micropoint, double-armed	994996
	5 - 0	1	75 cm	350 Mi.	13 mm micropoint, double-armed	994095


Polypropylene blue



Needle code	USP	EP	CM	Product No. (REF)
1/2 CIRCLE - TAPER POINT (ROUND-BODIED)				
2 x HR-13 	6 - 0	0.7	75 cm	9861
	5 - 0	1	75 cm	9862
	5 - 0	1	60 cm	9862S
2 x HR-17 	5 - 0	1	90 cm	9872N
	4 - 0	1.5	90 cm	9873N
	3 - 0	2	90 cm	9874N
2 x HR-20 	4 - 0	1.5	90 cm	98753N
	3 - 0	2	90 cm	98754N
2 x HR-22 	4 - 0	1.5	90 cm	9883N
2 x HR-25 	4 - 0	1.5	90 cm	98853N
	3 - 0	2	90 cm	98854N
	2 - 0	3	90 cm	98855N
3/8 CIRCLE - TAPER POINT (ROUND-BODIED)				
2 x TR-13 	6 - 0	0.7	75 cm	9991
	5 - 0	1	75 cm	9992
	5 - 0	1	90 cm	9992N
1/2 CIRCLE - CUTTING TAPER				
2xHSR 18 	4 - 0	1.5	90 cm	9943N
	3 - 0	2	90 cm	9944N
	2 - 0	3	90 cm	9945N
2xHSR 25 	4 - 0	1.5	90 cm	99483N
	3 - 0	2	90 cm	99484N
	2 - 0	3	90 cm	99485N

Surgical Stainless Steel (Sternotomy)

A special stainless steel monofilament made for the approximation of the sternum, very smooth and flexible with high tensile strength. The sterile Assut sternum wire is available with HCS 50 or TCS 50 needle on diameters USP -5-, -6- & -7-. Most modern packaging system for easy OR-handling.



Surgical Stainless steel monofilament

Needle code	USP	EP	Multipack	Product No. (REF)
1 x HCS-50	-7-	9	4 x 45 cm	CV 7605x4
	-6-	8	4 x 45 cm	CV 7624x4
	-6-	8	6 x 45 cm	CV 7622x4
	-6-	8	4 x 75 cm	CV 7623x4
	-5-	7	4 x 45 cm	CV 7601x4
	-5-	7	6 x 45 cm	CV 7620x4
	-5-	7	4 x 75 cm	CV 7621x4
	-5-	7	2 x 75 cm	CV 7802x2
1xTCS-50	-7-	9	4 x 45 cm	CV 7506x4
	-6-	8	4 x 45 cm	CV 7503x4
	-6-	8	6 x 45 cm	CV 7504x4
	-6-	8	4 x 75 cm	CV 7507x4
	-5-	7	4 x 45 cm	CV 7501x4
	-5-	7	6 x 45 cm	CV 7502x4
	-5-	7	4 x 75 cm	CV 7508x4



Temporary Cardiac Pacing Wire

A very fine braided multifilament stainless steel wire of excellent flexibility with coloured Polyethylen insulation, easy to identify. A special straight cutting needle RSK 90 or RSK 58 with breakaway feature, combined with a round - bodied needle (25 or 17 mm) or with a cutting taper needle (18 or 25 mm) or loop.

Needle code	USP	EP	CM	Colour	Product No. (REF)
HR 25 + RSK 90	2-0	3	60	black	CVP 150
	2-0	3	60	blue	CVP 151
HR 17 + RSK 90	3-0	2	60	blue	CVP 660
HR 17 + RSK 58	3-0	2	60	blue	CVP 656
	2-0	3	60	blue	CVP 658
HSR 18 + RSK 58	3-0	2	60	blue	CVP 652
HSR 25 + RSK 58	3-0	2	60	blue	CVP 653
LOOP + RSK 90	2-0	3	60	white	CVP 220
	2-0	3	60	blue	CVP 221



Assut Microsurgical Knives

Made from first quality stainless steel, the sterile Assut Microsurgical Knives ensure highest precision and exceptional sharpness. For further information refer to our special documentation.



Assut Knives

of Switzerland

Assut sutures

Assut Medical Sàrl
P.O. Box 5
CH - 1009 Pully-Lausanne
Switzerland
Tel. ++41-21-728.64.35
Fax. ++41-21-728.64.49
Email : info@assutsutures.com
www.assutsutures.com

Assut Knives