


## Single use Electrosurgical Pencil



### C-PEN

The CIMPAX<sup>®</sup> Single Use Electrosurgical Pencil C-PEN is a high quality pencil. The electrode is securely fitted with a hexagonal connection securing the electrode in the required position. C-PEN can be used with any electrosurgical generator available.

- Comfortable design and light weight
- Colour marked buttons
- Hexagonal electrode connection
- 3 meter cable, 2.35mm electrode
- 3 pin US-connector
- Single use
- Packaging: Box of 25 pcs, Sterile

Ref. code C31-08-001

CE0543