






## Post Operative Pain Management Simple and effective




CIMPAX® Bolus and Infiltration Kit  
Complete local anaesthesia kit for post-surgical pain relief  
Including CIMPAX® C-CAT catheter


CIMPAX® Bolus and Infiltrations Kit C-BIK, is the complete kit for post-operative pain management. C-BIK contains syringes for infiltration of the surgical site during the procedure and the C-CAT catheter for distribution of anaesthetics to the surgical site post-operatively. The catheter can be used with infusion pumps as well as for bolus injections. For bolus injections the C-CAT catheter offers minimal resistance due to its unique design.




- Precise, direct and concentrated pain relief
- Reduces or removes the use of narcotics
- Less or no side effects caused by narcotics
- Less pain (improved VAS score)
- Easier patient mobilisation
- Faster recovery, shortens hospital stay
- Increased patient satisfaction




Infiltration of the surgical site



Introduction of the catheter sheath



Introducing catheter through the sheath



Closing the wound and fixation of catheter

## Ordering information

### C-CAT

Catheter with 6 holes covering 2,5 cm	Ref. C44-21-005
Catheter with 12 holes covering 6 cm	Ref. C44-21-010
Catheter with 30 holes covering 15 cm	Ref. C44-21-015
Catheter with 60 holes covering 30 cm	Ref. C44-21-020
Packaging: Box of 14 kits, sterile	

### C-CAT, light

Catheter with 12 holes covering 6 cm	Ref. C44-21-110
Catheter with 30 holes covering 15 cm	Ref. C44-21-115
Packaging: Box of 14 kits, sterile	

### C-BIK

Complete Bolus and Infiltration Kit	
Catheter with 12 holes covering 6 cm	Ref. C44-21-001
Packaging: Box of 7 kits, sterile	