





C-VAC TORNADO™ Smoke Evacuation Electrosurgical Pencil



CIMPAX® C-VAC TORNADO™
Electrosurgical Pencil with Smoke Evacuation
LATEX, LEAD, PVC and PHTHALATE FREE

CIMPAX[®] C-VAC TORNADO[™] the ultimate combination of electrosurgery and removal of hazardous smoke in one instrument. The pencil has a unique triangular design which combined with the silicone cover offers the surgeon perfect ergonomic and safety grip. The ultra soft micro switch push buttons gives further comfort for the user due the exceptional tactile feedback. Since the C-VAC TORNADO[™] is Latex, lead, PVC and phthalate free you can perform surgery without concerns for your health, your colleagues and patient's health or the environment.

- Unique triangular design from tip to end
- Ultra soft push buttons for enhanced comfort
- Maximum safety grip silicone cover
- Transparent tip for clear view of the electrode
- LATEX, LEAD, PVC and PHTHALATE FREE
- Flexible tubing with integrated ESU cable
- Ready to use, no assembling of suction tip or electrode required.



Ordering information

C-VAC TORNADO[™] Integrated smoke evacuation pencil
4 meter lead, PVC and phthalate free ESU cable with
US-3-pin connector and 3 meter PVC and phthalate free
tubing with combined 10/22mm connector. Single use.

With 70 mm standard blade electrode

Ref. C31-19-001

With 70 mm coated blade electrode

Ref. C31-19-101

Packaging: Box of 10 pcs., sterile.